




Sanitising Wipes

Cleanline® Sanitising Wipes are a fast-acting cleaning and disinfecting wipe effective against a range of micro-organisms such as bacteria and yeast.

- Designed to meet the latest regulatory requirements for Maximum Residual Levels (MRLs) for surface disinfectants in food preparation and handling environments.
- Will clean and disinfect surfaces without using quats, PHMB or alcohol.
- Proven to kill micro organisms such as bacteria & yeasts.
- HFA approved.
- Independently tested to a range of European standards according to EN14885:2015.

DIRECTIONS FOR USE:

Remove lid and thread the first wipe through hole in centre of the lid. Replace the lid and tear wipes at perforations. Close the centre cap after use. For heavily soiled surfaces, more than one wipe may be required. To disinfect allow required contact time.



MICROBIOLOGICAL ACTIVITY

Tested independently in accordance with;

EN 1276 Minimum contact time of 1 minute.	Bactericidal action against <i>Pseudomonas aeruginosa</i> , <i>Staphylococcus aureus</i> , <i>Enterococcus hirae</i> , <i>Escherichia coli</i> , <i>Lactobacillus brevis</i> , <i>Serratia marcescens</i> , <i>Shigella sonnei</i> , <i>Salmonella typhimurium</i> , <i>Listeria monocytogenes</i> , <i>E. coli</i> 0157 and <i>Campylobacter jejuni</i> .
EN 13727 Minimum contact time of 1 minute.	Bactericidal action against <i>Enterococcus hirae</i> , <i>Pseudomonas aeruginosa</i> and <i>Staphylococcus aureus</i> .
EN 13697 Minimum contact time of 1 minute.	Bactericidal action against <i>Pseudomonas aeruginosa</i> , <i>Staphylococcus aureus</i> , <i>Enterococcus hirae</i> and <i>Escherichia coli</i> .
EN 16615 Minimum contact time of 5 minutes.	Bactericidal action against <i>Pseudomonas aeruginosa</i> , <i>Staphylococcus aureus</i> and <i>Enterococcus hirae</i> .
EN 1650 Minimum contact time of 1 minute.	Yeasticidal action against <i>Candida albicans</i> .
EN 16615 Minimum contact time of 1 minute.	Yeasticidal action against <i>Candida albicans</i> .

A summary of the tests is available on request.



CODE AND CASE SIZE:

- CL4028** - 200 Wipe Tub
- CL4029** - 1000 Wipe Bucket Tub
- CL4030** - 1500 Wipe Bucket Tub

QUALITY ASSURANCE:

This product is manufactured in the UK for Bunzl UK & Ireland and produced under ISO 9001 Quality Management System & ISO 14001 Environmental Management Systems.

ENVIRONMENTAL:

All surfactants used in Cleanline products comply with the current European regulations concerning biodegradability & protection of the environment.

INGREDIENTS & FORMULATION:

Specially formulated solution containing N,N-bis (3-aminopropyl) dodecylamine. The disinfectant is characterised by good activity against mycobacteria, low efficiency loss under protein load and little foaming as well as being quat, alcohol and PHMB free which makes it ideal for use within the food industry.

SAFE HANDLING:

Avoid contact with eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Dispose of wipes in accordance with local legislation. Do not flush wipes. Keep out of reach of children.

SAFE STORAGE:

Store sealed in original container in a cool, well ventilated area away from direct sunlight and extremes of temperature.

PO BOX 15247 Birmingham B23 3HN
Email: info@prime-source.co.uk

ISSUE DATE: 04/02/2020



Cleanline® is an exclusive brand of Bunzl UK & Ireland

© 2020 Bunzl UK Ltd. All Rights Reserved.